



FRIOS® | **DENTSPLY**
FRIADENT

Perfect light!
Perfect result!

FRIOS® straight and contra-angle handpieces with LED light ensure ideal intra-oral illumination thus creating the best conditions for perfect treatment results.

The most important advantages at a glance:

- Operate in daylight quality: Neutral white-light with high luminosity and a wider field of illumination than conventional halogen lamps.
- Independent electric power supply: Electricity for the LED light comes from the independent generator integrated in the instrument. Irrespective of the drive system you work with.
- Compatible with every motor with ISO coupling
- Extremely low heat generation
- Longer service life than conventional halogen lamps
- Sterilizable to 135 °C
- Thermal disinfectable



FRIOS® MicroSaw acc. to Prof. Dr. F. Khoury

NEW



FRIOS® straight and contra-angle handpieces with LED light

DENTSPLY Friadent presents a team perfectly coordinated in function and handling – the FRIOS® Unit S/i and the new FRIOS® LED handpieces and contra-angle handpieces specially adapted to the hand motor, simple to use, reliable in operation and convincing with their results.

Technical data:

	FRIOS® WI-75 LED G Contra-angle Handpiece	FRIOS® SI-11 LED G Straight Handpiece
Order number	90-6810	90-6850
Translation	20 : 1	1 : 1
Coupling system	ISO 3964 DIN 13.940	ISO 3964 DIN 13.940
Light system	LED	LED
Power supply	via independent generator (G)	via independent generator (G)
Luminosity	e.g. 25,000 Lux* at 800 rpm	e.g. 25,000 Lux* at 16.000 rpm
Light color	Daylight quality	Daylight quality
Sterilizable	Autoclave, to 135 °C	Autoclave, to 135 °C

* directly on LED

Order no. 6-302031/001

Friadent GmbH · Postbox 71 01 11 · 68221 Mannheim/Germany
Phone: +49 621 4302-1265 · Fax: +49 621 47 13 06 · e-mail: info@friadent.de

Manufacturer FRIOS® Handpieces, FRIOS® Unit S/i:
W&H Dentalwerk Bürmoos GmbH · 5111 Bürmoos/Austria

CE 0297

www.dentsply-friadent.com

DENTSPLY
FRIADENT